



Aesthetics - MUSTER

SLIMPRESS

professional equipment for sequential press therapy treatments



muster

✉ bbazar@globelife.com

☎ +39 0331 1706328



SLIMPRESS

professional equipment for sequential press therapy treatments

Exerts on the lymphatic circuit a beneficial toning and draining that enhances the natural formation of nutrients, improving both the tone that the appearance of the skin. The natural increase in diuresis, patients with fluid retention, produces beneficial effects in support of elimination of waste retained in the tissues.

Leg 7 full elements bracelets:

- 2 inflatable units foot/ankle
 - 2 inflatable units legs
 - 1 inflatable unit abdomen
 - 2 inflatable units arms
- Quick closure velcro high resistance.

MAIN FEATURES

- Large color Touch Screen Display
- Compact and silent
- Intuitive and easy to use

4 pre-set operational programmes:

A-sequential from foot to arms, with inflation of individual sectors. Lightweight, ideal for maintenance.

B-intensive sequential from foot to arms, with inflation a couple of areas. For evident presence of cellulite, fat and generic water retention.

C-sequential from foot to arms, with maintenance of areas inflated. Intensive, concentrated on the abdomen and thighs.

D-sequential operating arms and legs.